



Instrument tables by Rodenstock

**EPN, EPN V**

**Single tables - the practical integration solution.  
Instrument table EPN.**

With an EPN-table you possess a very functional instrument carrier of excellent quality - the ideal instrument table for your Rodenstock Instruments. The tables require only minimum space and are made with a silently operating motorised height adjustment feature. As the elevating column is made of aluminium the EPN instrument carrier has a minimum weight. The base plate has an abrasion-resistant powder coating making the surface damage resistant and keeps this "new" look. The carriers of the EPN series are ergonomically designed and have modern look. Depending on your needs two types of scratch resistant table tops are available: EPN instrument table with a rectangular table top for one instrument and the EPN-V instrument carrier for two instruments. The tables are available with two types of table bases, casters for easy moving of the unit or with a flat base for a fixed location.



**Technical specifications:**

<b>Base Module</b>	
Adjustable Height	290 mm
Height including desktop	min. 570 mm, max. 860 mm
Weight	32 kg
Nominal voltage	230 V
Power Consumption	210 W

<b>Modell</b>	<b>EPN</b>	<b>EPN V</b>
Length	500 mm	560 mm
Width	600 mm	1000 mm
Weight	7 kg	9 kg

**Ordering information.**

Please order one base module in conjunction with one desk top.

Subject to changes due to technical development.

For further information on EPN instrument carrier and Rodenstock Instruments please call us.

We would be pleased to inform you.

**WECO Optik GmbH**  
Administration Jägerstraße 58  
D-40231 Düsseldorf  
Telefon +49-211-21 04-104  
Telefax +49-211-21 04-251  
info@weco-instruments.com  
www.weco-instruments.com

Distributed by: